

# SPSeq-Vazagabiome 1

## Materials and Methods (continue)

Bacterial strains used in this study (Ref strains) add other ID#

- ATCC 8014 *Lactobacillus plantarum*
- ATCC 6633 *Bacillus subtilis*
- DSM 17938 *Lactobacillus reuteri*
- DHS2 Inulinogen saccharifying efficiency comp. cells Cat# 18265017

Note: A 25% Glycerol stock was prepared, for all the bacterial strains listed above.  
 New box of bacterial strains Pelekanos' lab has established  
 Labeling was done the following BVP# 1 b = bacteria  
 2 VP = Vazagabiome Pelekanos

## Growth conditions *Lactobacillus* sp (Gram +)

*Lactobacillus* strains ATCC 8014 and DSM 17938 were grown @ 37°C in MRS Broth. O/N cultures (1 volume) not exceeding OD 1 were used to start 100ml (MRS Broth) cultures. OD<sub>600</sub> 0.05 - 0.07

Bacteria was grown to log / mid-log phase ~ 5h for ATCC 8014  
 ~ 7h for DSM 17938

OD<sub>600</sub> for biological replicates

	(OD <sub>600</sub> )		(OD <sub>600</sub> )
ATCC 8014 - 1	0.49	DSM 17938 - 1	0.34
ATCC 8014 - 2	0.47	DSM 17938 - 2	0.56
ATCC 8014 - 3	0.55	DSM 17938 - 3	0.376

Note: DSM 17938 needed ~ 7h to reach OD<sub>600</sub> log-phase  
 ATCC 8014 needed ~ 5h to reach OD<sub>600</sub> log-phase

## Growth conditions *Bacillus subtilis* and DHS2 (Gram + Control Strains)

*Bacillus subtilis* ATCC 6633 and DHS2 were grown in LB @ 37°C  
 O/N cultures were used to start 60ml (LB) cultures to OD<sub>600</sub> 0.05 - 0.07  
 Bacteria were grown to mid-log phase ~ 3h

OD<sub>600</sub> for three biological replicates of control strains

ATCC 6633-1	<sup>OD<sub>600</sub></sup> 0.8	DH52-1	<sup>OD<sub>600</sub></sup> 0.88
ATCC 6633-2	0.78	DH52-2	0.89
ATCC 6633-3	0.98	DH52-3	0.7

Note ATCC 6633 B Substratum exhibited flocculent growth.  
possibly because it was grown @ 37°C. (NC @ 30°C)

Both strains needed ~ 3h to reach mid-log phase

### Culture Conditions Continue

20ml of Culture was harvested for each treatment  
and pellet was frozen in EtOH dry ice bath for  
RNA analysis

Experiments were done in three biological replicates  
(see <sup>OD<sub>600</sub></sup> page 4, 5)

For strains ATCC 8014 and DSM 17938 I had total of  
100ml for each replicate. (in FKS Broth) 5 conditions

2x Control (untreated) 1x CATI 1x SHX and 1x HUP

For Control strains ATCC 6633 and DH52 I had total  
of 60ml culture for each replicate in LB 3 conditions

2x Control (untreated) 1x CATI

### Chloramphenicol (CATI) treatments

CATI was added to final concentration of 100 µg/ml

(20ml Culture + 67 µl of 30mg/ml stock in EtOH)

incubated for 5' @ 37°C shaker and harvested

on ice (50ml conical filled with ice) containing

67 µl CATI in table top centrifuge for 2' full speed  
Pellet was shock frozen and stored for RNA analysis

## DL-Serine hydroxamate treatments (SHX)

SHX was added to final concentration of 1.25 mg/ml  
(20ml culture + 500  $\mu$ l of 50 mg/ml SHX stock in H<sub>2</sub>O)

10' @ 37°C, harvested in table-top centrifuge full speed  
2'. Pellet was shock frozen and stored @ -20°C for  
RNA analysis

## Mupirocin (MUP) treatments

MUP was added to final concentration of 65  $\mu$ g/ml  
(130  $\mu$ l of MUP 10 mg/ml stock in DMSO to 20ml culture)

## Stationary grown bacteria cultures

for stationary phase cultures strain ATCC 8014 was  
grown in MRS Broth for 27h @ 37°C

OD<sub>600</sub> for 3 biological replicates:

ATCC 8014-S1	4.3
ATCC 8014-S2	4.5
ATCC 8014-S3	4.65

2ml culture was harvested and stored for RNA analysis

Note: prior I setup cultures for stationary phase  
and grew them 72h for strain ATCC 8014 and DSM 17938.  
RNA of these samples was degraded. Therefore I decided  
to only use ATCC 8014 and left out DSM 17938 due to  
slow growth.

## Heat shock and nutrient deprivation culture conditions

*Lactobacillus plantarum* ATCC 8014 was grown in MRS Broth

@ 37°C. I used the 27h grown stationary phase

Cultures as precultures (see above) and inoculated in

40ml MRS Broth to OD<sub>600</sub> 0.05. Bacteria was grown to

log-phase (OD<sub>600</sub>) ~ 4.5h

ATCC 8014-4	OD <sub>600</sub> 0.2
ATCC 8014-5	0.31
ATCC 8014-6	0.39

10ml of each replicate for untreated control was harvested and Pellet was stored for RNA analysis.

45ml of each replicate was harvested after corresponding stress condition

### Heat Shock treatments HS

15ml Culture was harvested in table-top centrifuge

@ full speed for 2'

ATCC 8014-HS1 (4)

ATCC 8014-HS2 (5)

ATCC 8014-HS3 (6)

Note: (1ml <sup>HB</sup> Broth was preheated in Thermomixer to 60°C / in 2ml Epi tube)

Pellet was resuspended in preheated media and cells were incubated for 15' @ 60°C light shaking in Thermo Mixer

Pellet was stored for RNA analysis

### Nutrient Deprivation treatments (ND) 45' in 0.5xLB

15ml Culture was harvested in table-top centrifuge full speed

2'. Pellet was washed extensively (50x) 0.5xLB

→ Start timer

Spin 2' full speed. Sup was ~~disregarded~~ <sup>discarded</sup> pellet was resuspended in 20ml 0.5xLB media

Culture was incubated @ 37°C shaking for remaining time (45' total start of timer @ wash)

Culture was harvested and Pellet stored for RNA analysis.

ATCC 8014-ND1 (4) & # corresponds to to <sup>biological</sup> replicate control sample

ATCC 8014-ND2 (5)

ATCC 8014-ND3 (6)

Note: I called it Nutrient deprivation, however I am not sure this is valid because LB is rich media. Talked to recent he recommends to call it unfavorable culture media condition

## RNA Isolation (Coffey Protocol Book page 63 glass bead RNA iso1)

RNA was isolated as described in Coffey Protocol book (see attached PDF) with minor modifications)

- Step 4 vortex in multimixer for 90"
- Step 6 vortex for 90" and spin for 2' full speed
- Step 7 vortex for 30" spin for 1' full speed
- Step 8 vortex for 30" spin for 1' full speed
- Step 10 centrifuge for 20' full speed 4°C
- Step 11 resuspend pellet/RNA in 27µl of nuclease-free H<sub>2</sub>O

Quantification of RNA was done by Nanodrop  
 ↳ see attachments, Vagagabiome 1, RNA quantification

used the 1-10  
 dilutions to determine  
 final conc



5P Seq (20x11)

• A 1µg 1ml working RNA stock was prepared for library prep. of all the RNA samples. Original RNA stock is stored @ -80

• 1µg of total RNA was run on 1.2% Agarose gel in TAE to assess quality of RNA preps. total of 60 RNA preps

1A TOP RNApreps 1-1 - 1-16  
 User 2018-05-21\_18h23m23s

2A TOP RNApreps 3-1 - 3-16  
 User 2018-05-23\_19h55m32s



1B bottom RNApreps 2-1 - 2-16

2B bottom RNApreps 4-1 - 4-12

Note: 1µg of total RNA separates roughly to: 600ng 23S  
 300ng 16S  
 100ng rRNA

Note: Quality of RNA was OK

- little smear/degradation of sample 1-5
- for samples 4-3, 1-14, 1-16 23S/16S Ratio not great
- samples 4-4, 4-5, 4-6 low yield of RNA (see Nanodrop concentrations PDF)



Legend for RNA quality check for biological replicates #1, #2, #3, #4

1A	E#1	1B	E#2	2A	E#3
1-1	ATCC 8014-1 Chr	2-1	ATCC 8014-2 Chr	3-1	ATCC 8014-3 Chr
(1-2)	ATCC 8014-1 Chr	(2-2)	ATCC 8014-2 Chr	(3-2)	ATCC 8014-3 Chr
1-3	ATCC 8014-1 CAM	2-3	ATCC 8014-2 CAM	3-3	ATCC 8014-3 CAM
1-4	ATCC 8014-1 SHX	2-4	ATCC 8014-2 SHX	3-4	ATCC 8014-3 SHX
1-5	ATCC 8014-1 MUP	2-5	ATCC 8014-2 MUP	3-5	ATCC 8014-3 MUP
(1-6)	DSM 17938-1 Chr	2-6	DSM 17938-2 Chr	3-6	DSM 17938-3 Chr
(1-7)	DSM 17938-1 CH	(2-7)	DSM 17938-2 CH	(3-7)	DSM 17938-3 CH
1-8	DSM 17938-1 CAM	2-8	DSM 17938-2 CAM	3-8	DSM 17938-3 CAM
1-9	DSM 17938-1 SHX	2-9	DSM 17938-2 SHX	3-9	DSM 17938-3 SHX
1-10	DSM 17938-1 MUP	2-10	DSM 17938-2 MUP	3-10	DSM 17938-3 MUP
1-11	ATCC 6633-1 Chr	2-11	ATCC 6633-2 Chr	3-11	ATCC 6633-3 Chr
(1-12)	ATCC 6633-1 CH	(2-12)	ATCC 6633-2 CH	(3-12)	ATCC 6633-3 CH
1-13	ATCC 6633-1 CAM	2-13	ATCC 6633-2 CAM	3-13	ATCC 6633-3 CAM
1-14	DH52-1 Chr	2-14	DH52-2 Chr	3-14	DH52-3 Chr
(1-15)	DH52-1 CH	(2-15)	DH52-2 CH	(3-15)	DH52-3 CH
1-16	DH52-1 CAM	2-16	DH52-2 CAM	3-16	DH52-3 CAM

## 2B E#4 Station samples

- 4-1 ATCC 8014-4 Chr
- 4-2 ATCC 8014-5 Chr
- 4-3 ATCC 8014-6 Chr
- 4-4 ATCC 8014-HS1 (4)
- 4-5 ATCC 8014-HS2 (5)
- 4-6 ATCC 8014-HS3 (6)
- 4-7 ATCC 8014-ND1 (4)
- 4-8 ATCC 8014-ND2 (5)
- 4-9 ATCC 8014-ND3 (6)
- 4-10 ATCC 8014-S1 (4)
- 4-11 ATCC 8014-S2 (5)
- 4-12 ATCC 8014-S3 (6)

NOTE:

- Chr = Control (untreated)
- CAM = Chloramphenicol 100 µg/mL
- SHX = DL Serine hydroxamate 125 µg/mL
- MUP = Mupirocin 60 µg/mL
- HS = HEAT SHOCK 15' @ 60°C
- ND = Nutrient deprivation 15' @ 0.3% OD ~ 4
- S = Stationary phase 27h

Additional Libraries to prepare for MIXED POPULATION Tgs

- 5-1 = 2-6 + 2-1 1:1
- 5-2 = 2-6 + 2-3 1:1
- 5-3 = 2-6 + 2-4 1:1
- 5-4 = 2-6 + 2-5 1:1
- 5-5 = 2-6 + 2-1 1:10
- 5-6 = 2-6 + 2-3 1:10
- 5-7 = 2-6 + 2-4 1:10
- 5-8 = 2-6 + 2-5 1:10
- 5-9 = 2-6 + 2-1 1:100
- 5-10 = 2-6 + 2-3 1:100
- 5-11 = 2-6 + 2-4 1:100
- 5-12 = 2-6 + 2-5 1:100

5 µl + 5 µl

9 µl + 1 µl

9.9 + 0.1 µl

- See modified SPSeq Protocol Ins PDF

NOTE: total of 72 Libraries; 60 x Standard SPSeq  
12 x random frag control

### Random frag controlsamples

E1 1-2, 1-7, 1-12, 1-15  
E2 2-2, 2-7, 2-12, 2-15  
E3 3-2, 3-7, 3-12, 3-15

Onase treatment as described in Protocol.

First round was E1, E2, E3, E4 (except random frag control samples)

Left the samples on @ -20

next day I did E5 and the random frag controls (Pink etc)

RNA Ligation was done for all of them on @ 16C (17h)

next day, Etho precipitation and elute in 16p nucleic acid

### Rabbit RNA

- for RNA depletion RNA concentration cannot exceed 1pg total. Unf. We do not have enough RNA rabbit reagent

Therefore I can only measure a couple and decide then what I do.

### Rabbit measurements

sample	ng/ml	ng/ml	total in 16p
1-1	0.186	140	1971
1-6	0.06	60	800
1-11	0.19	190	2565
1-14	0.07	70	945
2-1	0.1	100	1350
2-6	0.07	70	945
2-11	0.19	190	2565
2-14	0.06	60	810
3-1	0.07	70	945
3-6	0.17	170	2295
3-11	0.03	30	405
3-14	0.03	30	405
4-1	0.17	170	2295
4-2	0.12	120	1620
4-5	0.18	180	2430

sample	ng/ml	ng/ml	total in 16p
5-1	0.19	190	2565
5-2	0.19	190	2565
5-3	0.18	180	2430

I did compare ~~some~~ <sup>measure</sup> some samples with nano drop

but the values were all over the place

Unf. I cannot do it like that

	nano drop	Rabbit
2-9 1-1	277 ng/ml	146
1-6	180 ng/ml	60
1-11	400 ng/ml	190

Given that the Nanodrop values and Qubit values  
do not correlate and the Ribozero Kit guide specifically  
says that RNA concentrations should be measured with Qubit  
we have a problem: We do not have enough RNA Qubit  
solution.

I decided that I will take  $\frac{1}{2}$  V for rRNA depletion  
as it is stated in my protocol. For the reason that my  
highest concentration is 2.5  $\mu$ g and I want 400 ng in 13.5  $\mu$ l  
The rRNA depletion kit can go as low as 100 ng.

Start with E1 and E2 total of 32 samples

For bead washing I need 45  $\mu$ l/sample  $\times$  33 1.5 mL beads  
after wash <sup>batch</sup> resuspend in 400  $\mu$ l resusp. Buff + 16.5  $\mu$ l Rib. 6.

Do in 2 batches 250  $\mu$ l beads, wash, wash

resuspend in 200  $\mu$ l resp. Buff + 8.25  $\mu$ l pull them and add 16.5

mix 6  $\mu$ l nuclear free luo

+ 7.5 $\mu$ l RNA	}	6 $\mu$ l	148.5
+ 4.5 removal		MM for 33.5	66
+ 2 $\mu$ l 10x Ract. Buff			

Next set E3, E4, E5 total 40 / 41 sam

7.5 $\mu$ l	
4.5 rem	187
2 $\mu$ l 10x	83

batch wash strip beads	45 $\times$ 41	1845
		922.5 each

492 susp  
246 each  
20  $\mu$ l Rib. 6



3

Pls RT I will do Ampure bead purification and check

70 ml Sample

Topic of Physics HW

20 ml - 5 ml A

023 - 7

• binding to sheep beads.	M280	hydrophobic
	M270	hydrophobic

o O/n Com on adaptor (yehon stat 0 895)

June took my shift out @ 10am (213k total t. of L)

Ask  $\rightarrow$  Wash with ligden and, and let the beads settle as written

→ Wash with hydrogen and 1 chd bet the beads. Since  
 # For the PCR Amplification 1 chd total of 50pc  
 should be beads & if that fails 40 pc

• stored the beads with 100% EtOH in 4°C 22.06.18

FLR @ -20 will do Airport Preparation on Monday 25.

# Index and Multiplexbarcode for Library <sup>PCR</sup> amplification

14

	NEBL	PE2-MPX		NEBL	PE2-MPX
1-1	503	04	2-1	504	04
1-2	503	05	2-2	504	05
1-3	503	06	2-3	504	06
1-4	503	07	2-4	504	07
1-5	503	08	2-5	504	08
1-6	503	09	2-6	504	09
1-7	503	10	2-7	504	10
1-8	503	11	2-8	504	11
1-9	503	12	2-9	504	12
1-10	503	13	2-10	504	13
1-11	503	14	2-11	504	14
1-12	503	15	2-12	504	15
1-13	503	16	2-13	504	16
1-14	503	17	2-14	504	17
1-15	503	18	2-15	504	18
1-16	503	19	2-16	504	19

	NEBL	PE2-MPX		NEBL	PE2-MPX
3-1	505	04	4-1	506	04
3-2	505	05	4-2	506	05
3-3	505	06	4-3	506	06
3-4	505	07	4-4	506	07
3-5	505	08	4-5	506	08
3-6	505	09	4-6	506	09
3-7	505	10	4-7	506	10
3-8	505	11	4-8	506	11
3-9	505	12	4-9	506	12
3-10	505	13	4-10	506	13
3-11	505	14	4-11	506	14
3-12	505	15	4-12	506	15
3-13	505	16			
3-14	505	17			
3-15	505	18			
3-16	505	19			

	NEBL	PE2-MPX
5-1	503	20
5-2	503	21
5-3	503	22
5-4	504	20
5-5	504	21
5-6	504	22
5-7	505	20
5-8	505	21
5-9	505	22
5-10	506	20
5-11	506	21
5-12	506	22